

Using AED for Adults

<p>1. Turn on AED and follow the voice prompts.</p>	<p>5. Clear for analysis.</p> 
<p>2. Remove all clothing covering the chest, if necessary.</p>	
<p>3. Attach pads correctly.</p> <ul style="list-style-type: none">• Place one pad on upper right side of chest.• Place one pad on lower left side of chest, a few inches below the left armpit.• Pads should not touch.  	<p>6. Clear for shock.</p>
<p>4. Plug the pad connector cable into the AED, if necessary.</p> 	<p>7. Push shock button to deliver shock.</p>  <p>8. Immediately start CPR.</p> <p>Note: Skill practice ends here.</p> 